



www.krnano.com


안과

OPHTHALMIC INSTRUMENTS

Page 098 - 115




Weiss
06-001-01 Solid for adults
06-001-02 Solid for infants
06-001-03 Fenestrated for adults
06-001-04 Fenestrated for infants



Castroviejo
06-002-01 Small
06-002-02 Medium
06-002-03 Large




Holtzheimer 4-teeth
06-003-10 10 cm 4ⁿ Sharp Str.
06-004-10 10 cm 4ⁿ Blunt Str.
06-005-10 10 cm 4ⁿ Sharp Cvd
06-006-10 10 cm 4ⁿ Blunt Cvd




Barraquer
06-007-01 Fenestrated
06-007-02 Fenestrated for infants




Desmarres
9 cm / 3½




20mm
1/1
06-008-09



25mm
06-009-09



30mm
06-010-09



Lambert
9.5 cm / 3½
06-011-09



1/2

Bishop-Harmon
06-012-09 - 06-015-09
8.5cm / 3 1/2"



1/1

06-012-09 Fine Tip
06-013-09 Standard



1/1

1x2 teeth
06-014-09 Fine Tip
06-015-09 Standard



Gill
06-016-07 - 06-018-07
7.0cm / 3 3/4"



1/1

1x2 teeth
06-016-07 Extra Delicate
06-017-07 Standard



1/1

Gill-Chandler
06-018-07



1/2

Gill-Hess
06-019-07
7.5 cm / 3"




1/1


1x2 teeth




Colibri
With tying platform
8.3 cm / 3 1/4"




06-020-08



06-021-08
Extra delicate





Bonn
Suture Forceps
8.3 cm / 3 1/4"





06-022-08





Castroviejo
06-023-11
11cm / 4½"





Castroviejo
06-024-12 0.12 mm
06-024-30 0.30 mm
06-024-50 0.50 mm
06-024-70 0.70 mm
11 cm / 4½"




Castroviejo TC
06-025-30 0.30 mm
06-025-50 0.50 mm
10.5 cm / 4½"

Tip 0.3 mm



1/1
06-026-11




06-027-11

Berke
10.8 cm / 4¼"



Micro Suture Tying


½



1/1

06-028-12 12 cm / 4¾"
06-028-15 15 cm / 6"

Tip 0.6 cm




1/1

06-028-18 18 cm / 7"
06-028-23 23 cm / 9"



1/1

06-029-12 12 cm / 4¾"
06-029-15 15 cm / 6"



1/1

06-029-18 18 cm / 7"
06-029-23 23 cm / 9"



Rhoton
18 cm / 7"

½




1/1


06-030-18 Str
06-031-18 Cvd




1/1




06-032-11
11 cm / 4½"




06-033-12
12 cm / 4¾"




06-034-09
09 cm / 3½"




06-035-09
09 cm / 3½"




Westcott
06-036-11 - 06-037-11
11 cm / 4½"



06-037-11



06-038-18
06-040-18
18.5 cm / 7½"



06-038-18



06-039-18



06-040-18



Vannas-Tubingen
06-041-08 Str
06-042-08 Cvd
08 cm / 3"



06-043-08 Str



06-044-08 Cvd



06-045-08 Bent



06-046-08 Angled on flat




06-047-08 Str



06-048-08 Cvd

Vannas very delicate
08 cm / 3"

Traberular
08 cm / 3"




06-049-18

06-050-18



1/2

06-049-18 - 06-050-18
18 cm / 7 1/4"



06-051-18



1/2


Hepp/Scheidel
18 cm / 7 1/4"



06-052-14
14 cm / 5 1/2"



06-053-15
15 cm / 6"




06-054-15
15 cm / 6"
1 blade serr.
Reill




1/2

06-052-14 - 06-054-15


Ribbon Scissors




06-055-10
10 cm / 4"




06-056-10
10 cm / 4"




06-057-05
9.5 cm / 3 3/4"




06-058-09
9.5 cm / 3 3/4"




Stevens
06-059-10
10 cm / 4"



Stevens
06-060-10
10 cm / 4"




Stevens
06-061-10
10 cm / 4"




Stevens
06-062-10
10 cm / 4"







2/3
06-063-09
9.5 cm / 3 3/4"




2/3
06-064-09
9.5 cm / 3 3/4"



2/3
06-065-09
9.5 cm / 3 3/4"




2/3
06-066-09
9.5 cm / 3 3/4"




2/3
Knapp
06-067-12
12 cm / 4 3/4"




2/3
Knapp
06-068-12
12 cm / 4 3/4"




2/3
Knapp
06-069-12
12 cm / 4 3/4"



2/3
Knapp
06-070-12
12 cm / 4 3/4"




06-071-10 TC
10 cm / 4"




06-072-10 TC
10 cm / 4"




06-073-09 TC
9.5 cm / 3¾"




06-074-09 TC
9.5 cm / 3¾"




Stevens TC
06-075-10
10 cm / 4"



Stevens TC
06-076-10
10 cm / 4"




Stevens TC
06-077-10
10 cm / 4"



Stevens TC
06-078-10
10 cm / 4"






1/1

2/3


06-079-09 TC
9.5 cm / 3¾"



1/1

2/3


06-080-09 TC
9.5 cm / 3¾"



1/1

2/3

06-081-09 TC
9.5 cm / 3¾"



1/1

2/3


06-082-09 TC
9.5 cm / 3¾"



1/1

2/3


Knapp TC
06-083-12
12 cm / 4¾"



1/1

2/3


Knapp TC
06-084-12
12 cm / 4¾"



1/1

2/3


Knapp TC
06-085-12
12 cm / 4¾"




1/1

2/3


Knapp TC
06-086-12
12 cm / 4¾"




06-087-10 SC
10 cm / 4"




06-088-10 SC
10 cm / 4"




06-089-09 SC
9.5 cm / 3 3/4"




06-090-09 SC
9.5 cm / 3 3/4"




Stevens
06-091-10 SC
10 cm / 4"



Stevens
06-092-10 SC
10 cm / 4"




Stevens
06-093-10 SC
10 cm / 4"



Stevens
06-094-10 SC
10 cm / 4"






1/1

2/3


06-095-09 SC
9.5 cm / 3¾"



1/1

2/3


06-096-09 SC
9.5 cm / 3¾"



1/1

2/3

06-097-09 SC
9.5 cm / 3¾"



1/1

2/3


06-098-09 SC
9.5 cm / 3¾"



1/1

2/3


Knapp SC
06-099-12
12 cm / 4¾"



1/1

2/3


Knapp SC
06-100-12
12 cm / 4¾"



1/1

2/3

Knapp SC
06-101-12
12 cm / 4¾"



1/1

2/3


Knapp SC
06-102-12
12 cm / 4¾"




Micro-Castroviejo
With out lock
9.5 cm / 3½"



Micro-Castroviejo
With lock
9.5 cm / 3½"




Barraquer
Without lock
12 cm / 4¾"



Micro-Castroviejo TC
With out lock
14 cm / 5½"



Micro-Castroviejo TC
With lock
14 cm / 5½"



Barraquer
With lock
12 cm / 4¾"





Castroviejo
Without lock
18 cm / 7"

Delicate 1/1



- 06-115-14 14 cm / 5½" Str/Smooth
- 06-115-18 18 cm / 7" Str/Smooth
- 06-115-21 21 cm / 8¼" Str/Smooth
- 06-115-23 23 cm / 8¾" Str/Smooth
- 06-116-14 14 cm / 5½" Str/Serrated
- 06-116-18 18 cm / 7" Str/Serrated
- 06-116-21 21 cm / 8¼" Str/Serrated
- 06-116-23 23 cm / 8¾" Str/Serrated

Micro 1/1



- 06-117-14 14 cm / 5½" Cvd/Smooth
- 06-117-18 18 cm / 7" Cvd/Smooth
- 06-117-21 21 cm / 8¼" Cvd/Smooth
- 06-117-23 23 cm / 8¾" Cvd/Smooth
- 06-118-14 14 cm / 5½" Cvd/Serrated
- 06-118-18 18 cm / 7" Cvd/Serrated
- 06-118-21 21 cm / 8¼" Cvd/Serrated
- 06-118-23 23 cm / 8¾" Cvd/Serrated



Castroviejo
With lock
18 cm / 7"

Delicate 1/1



- 06-119-14 14 cm / 5½" Str/Smooth
- 06-119-18 18 cm / 7" Str/Smooth
- 06-119-21 21 cm / 8¼" Str/Smooth
- 06-119-23 23 cm / 8¾" Str/Smooth
- 06-120-14 14 cm / 5½" Str/Serrated
- 06-120-18 18 cm / 7" Str/Serrated
- 06-120-21 21 cm / 8¼" Str/Serrated
- 06-120-23 23 cm / 8¾" Str/Serrated

Micro 1/1



- 06-121-14 14 cm / 5½" Cvd/Smooth
- 06-121-18 18 cm / 7" Cvd/Smooth
- 06-121-21 21 cm / 8¼" Cvd/Smooth
- 06-121-23 23 cm / 8¾" Cvd/Smooth
- 06-122-14 14 cm / 5½" Cvd/Serrated
- 06-122-18 18 cm / 7" Cvd/Serrated
- 06-122-21 21 cm / 8¼" Cvd/Serrated
- 06-122-23 23 cm / 8¾" Cvd/Serrated



Barraquer
Without lock
23 cm / 9"



1/1
Delicate



1/1
Micro

- 06-123-14 14 cm / 5½" Str
- 06-123-16 16 cm / 6" Str
- 06-123-18 18 cm / 7" Str
- 06-123-21 21 cm / 8¼" Str
- 06-123-23 23 cm / 8¾" Str
- 06-124-14 14 cm / 5½" Cvd
- 06-124-16 16 cm / 6" Cvd
- 06-124-18 18 cm / 7" Cvd
- 06-124-21 21 cm / 8¼" Cvd
- 06-124-23 23 cm / 8¾" Cvd

- 06-123-14 14 cm / 5½" Str
- 06-123-16 16 cm / 6" Str
- 06-123-18 18 cm / 7" Str
- 06-123-21 21 cm / 8¼" Str
- 06-123-23 23 cm / 8¾" Str
- 06-124-14 14 cm / 5½" Cvd
- 06-124-16 16 cm / 6" Cvd
- 06-124-18 18 cm / 7" Cvd
- 06-124-21 21 cm / 8¼" Cvd
- 06-124-23 23 cm / 8¾" Cvd



Barraquer
With lock
23 cm / 9"



1/1
Delicate



1/1
Micro

- 06-127-14 14 cm / 5½" Str
- 06-127-15 15 cm / 6" Str
- 06-127-18 18 cm / 7" Str
- 06-127-21 21 cm / 8¼" Str
- 06-127-23 23 cm / 8¾" Str
- 06-128-14 14 cm / 5½" Cvd
- 06-128-15 15 cm / 6" Cvd
- 06-128-18 18 cm / 7" Cvd
- 06-128-21 21 cm / 8¼" Cvd
- 06-128-23 23 cm / 8¾" Cvd

- 06-129-14 14 cm / 5½" Str
- 06-129-15 15 cm / 6" Str
- 06-129-18 18 cm / 7" Str
- 06-129-21 21 cm / 8¼" Str
- 06-129-23 23 cm / 8¾" Str
- 06-130-14 14 cm / 5½" Cvd
- 06-130-15 15 cm / 6" Cvd
- 06-130-18 18 cm / 7" Cvd
- 06-130-21 21 cm / 8¼" Cvd
- 06-130-23 23 cm / 8¾" Cvd






Castroviejo TC
Without lock
18 cm / 7"

- 06-131-14 14 cm / 5½" Str
- 06-131-18 18 cm / 7" Str
- 06-131-21 21 cm / 8¼" Str
- 06-131-23 23 cm / 8¾" Str
- 06-132-14 14 cm / 5¼" Cvd
- 06-132-18 18 cm / 7" Cvd
- 06-132-21 21 cm / 8¼" Cvd
- 06-132-23 23 cm / 8¾" Cvd



Castroviejo TC
With lock
18.0 cm / 7"

- 06-133-14 14 cm / 5½" Str
- 06-133-18 18 cm / 7" Str
- 06-133-21 21 cm / 8¼" Str
- 06-133-23 23 cm / 8¾" Str
- 06-134-14 14 cm / 5¼" Cvd
- 06-134-18 18 cm / 7" Cvd
- 06-134-21 21 cm / 8¼" Cvd
- 06-134-23 23 cm / 8¾" Cvd



Barraquer TC
With lock
23 cm / 9"



06-135-14	14 cm / 5½"	Str
06-135-15	15 cm / 6"	Str
06-135-18	18 cm / 7"	Str
06-135-21	21 cm / 8¾"	Str
06-135-23	23 cm / 8¾"	Str
06-136-14	14 cm / 5½"	Cvd
06-136-15	15 cm / 6"	Cvd
06-136-18	18 cm / 7"	Cvd
06-136-21	21 cm / 8¾"	Cvd
06-136-23	23 cm / 8¾"	Cvd

06-137-14	14 cm / 5½"	Str
06-137-15	15 cm / 6"	Str
06-137-18	18 cm / 7"	Str
06-137-21	21 cm / 8¾"	Str
06-137-23	23 cm / 8¾"	Str
06-138-14	14 cm / 5½"	Cvd
06-138-15	15 cm / 6"	Cvd
06-138-18	18 cm / 7"	Cvd
06-138-21	21 cm / 8¾"	Cvd
06-138-23	23 cm / 8¾"	Cvd




Barraquer TC
With lock
23 cm / 9"



06-139-14	14 cm / 5½"	Str
06-139-15	15 cm / 6"	Str
06-139-18	18 cm / 7"	Str
06-139-21	21 cm / 8¾"	Str
06-139-23	23 cm / 8¾"	Str
06-140-14	14 cm / 5½"	Cvd
06-140-15	15 cm / 6"	Cvd
06-140-18	18 cm / 7"	Cvd
06-140-21	21 cm / 8¾"	Cvd
06-140-23	23 cm / 8¾"	Cvd

06-141-14	14 cm / 5½"	Str
06-141-15	15 cm / 6"	Str
06-141-18	18 cm / 7"	Str
06-141-21	21 cm / 8¾"	Str
06-141-23	23 cm / 8¾"	Str
06-142-14	14 cm / 5½"	Cvd
06-142-15	15 cm / 6"	Cvd
06-142-18	18 cm / 7"	Cvd
06-142-21	21 cm / 8¾"	Cvd
06-142-23	23 cm / 8¾"	Cvd






Castroviejo
06-143-01 15 mm Str
06-143-02 15 mm Str with teflon
06-143-03 15 mm Cvd
06-143-04 15 mm Cvd with teflon




Castroviejo
06-144-01 20 cm Str
06-144-02 20 cm Str with teflon
06-144-03 20 cm Cvd
06-144-04 20 cm Cvd with teflon




06-145-08 8 mm
06-145-10 10 mm
06-145-12 12 mm
06-145-14 14 mm
06-145-16 16 mm
06-145-18 18 mm


Desmarres
Eye Retractor



ALM
4x4 prongs
7 cm / 2 3/4"



06-146-07 Sharp




06-147-07 Blunt




ALM
4x4 prongs
10 cm / 4"



06-148-10 Sharp



06-149-10 Blunt




Graefe
Iris Hook
12 cm / 4 3/4"



Bowman
Silver
13 cm / 5"

06-152-00 Set


Mgn		
06-152-01	000/000	0.4 / 0.6
06-152-02	00/0	0.7 / 0.8
06-152-03	1/2	0.9 / 1.0
06-152-04	3/4	1.3 / 1.4
06-152-05	5/6	1.5 / 1.6



Wilder
06-153-01 Fine
06-153-02 Medium
06-153-03 Strong




Graefe



06-154-12
12 cm / 4 3/4"



06-154-13
13 cm / 5"



Meyhoefer

06-155-00	0.5 mm
06-155-01	1.0 mm
06-155-02	1.5 mm
06-155-03	1.8 mm
06-155-04	2.0 mm
06-155-05	2.5 mm
06-155-06	3.0 mm
06-155-07	3.5 mm
06-155-08	4.0 mm
06-155-09	4.5 mm



초음파 세척기가 없는 세척방법 (손세척방법)

거품이 적게 발생하는 중성 pH(6-8.5)세제를 이용하여 부드러운 솔과 스폰지를 사용하여 손으로 의료기를 세척하세요.

렌즈가 달린 의료기구는. pH가 너무 알칼리성이거나 너무 산성이면 얼룩이나 구멍이 생길 수 있습니다.

의료기구는. 강모의 솔이나, 연마제를 사용하면 안 됩니다 약모의 브러쉬나 스폰지를 물에 잠긴 의료기구와 함께 사용해야 합니다.

의료기구는 증류수 또는 탈염수로 헹구어야 합니다.

그런 다음 보푸라기가 없는 천이나 뜨거운 공기로 건조시킵니다. 세척을 위해 식염수 대신 증류수를 사용하거나 헹굼 도구를 사용하면 의료기구의 수명을 수년 연장할 수 있습니다.

의료기구에 표시되어 있는 검은 글자는 레이저



빈 관형(루멘) 품목을 세척할 때는 가는 브러쉬를 이용하여 중성 pH(6-8.5)세제를 사용하여 부드러운 브러쉬로 세척합니다.